



COUNTERTOP LABORATORY FRIDGE-FREEZER COMBI SOLID DOOR 83.5L WHITE



- ✓ **Operating Temperature Range of +2°C to +10°C for Fridge and -6°C to -9°C for Freezer**
Conforms to MHRA guidelines
- ✓ **NEW Unique Anti Microbial Handle**
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Fitted Rollers**
For ease of movement
- ✓ **Wire Shelves and Basket**
Enables optimum air flow inside unit



The Pharmacy
Guild of Australia

SKU CODE:	PRODUCT CODE:	DESCRIPTION:
444442356	ISR27CAU	Countertop Laboratory Fridge-freezer Combi Solid Door 83.5L White



FULL SPECIFICATION – ISR27CAU COUNTERTOP LABORATORY FRIDGE-FREEZER COMBI SOLID DOOR 83.5L WHITE

Description	Countertop
SKU Code	444442356
EAN Code	5052263023567
GDPA Code	ISR27CAU
Gross Capacity (L)	83.5 (Fridge: 77, Freezer: 6.5)
Net Capacity (L)	80.5
Unit External Width (mm)	502
Unit External Depth (mm)	540
Unit External Height (mm)	660
Product Weight (kg)	25
Packed Width (mm)	620
Packed Depth (mm)	595
Packed Height (mm)	675
Packed Weight (kg)	27
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	415 x 370 x 540
External Finish	White
Internal Finish	White
Type of Plastic	SRL800 (HIPS/ABS Mix)
Temperature Range	Fridge: +2°C to +10, Freezer: -6°C to -9°C
Colour	White/Green
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	976
Replacement Door Seal	No
Open Door Sensor	No
Lock	Yes
Adjustable Feet	2 adjustable feet, 2 rear rollers
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	2
Number of Shelf Positions	Set
Shelf Size	1: 438 x 274 x 6, 1: 438 x 137 x 6
Shelf Material	Wire
Moulded Drawers	0
Adjustable Door Bins	0
Interior Light	No
Fan	No
Spark Free	Yes
Auto or Manual defrost	Fridge: Automatic, Freezer: Manual
Controller type	Digital, External
Climate Class	N
Refrigerant Type	R134a
CFC Free	Yes
HFC Free	No
HCFC Free	Yes
Australian Standard Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	230
Aux Output	Yes
Insulation Thickness (mm)	30
Insulation Type	Cyclopentane
Anti Microbial Handle	Yes
Tariff Code	8418219100