

LSFSF242BT-UK

FREESTANDING LABORATORY
FREEZER SOLID DOOR 242L WHITE



- ✓ **-18°C to -20°C Operating Temperature Range**
 Conforms to MHRA guidelines.
- ✓ **Lock Fitted as Standard**
 Maintains security and ensures the unit can only be opened by authorised personnel.
- ✓ **Wire Shelves**
 Enables optimum air flow inside the unit.
- ✓ **Unique Anti Microbial Handle**
 Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **3 Year Extended Warranty**
 When serviced by an approved service technician, the warranty is extended to 3 years.



Specifications

Capacity	242L Gross / 229L Net	Spark Free	Yes
Dimensions External	595w x 600d x 1560h	Internal Fan	No
Dimensions Internal	420w x 455d x 1345h	Temperature range	-18°C to -20°C
Lock	Yes, 2 Keys	Controller	Digital, External
Door	Solid	Interior	HIPS
Door Style	Right Hand Hinge	Housing	Powder coated sheet
Weight	58.2 kg	Voltage 50/60Hz	230V - 50Hz
Insulation Thickness	85mm	Energy Consumption	0.99kWh/24h
Sections	6	Internal light	No
Sections Material	HIPS Shelf with flap to open	Door Open Alarm	Yes
Defrost	Manual	Access Port	Yes - Rear Bottom Right (16mm)
Refrigerant	R600A / 64g	Noise Level (db)	36.5dB
Adjustable Feet (No & Position)	2 adjustable front feet	Anti Microbial Handle	Yes
Country of Origin	United Kingdom	Colour	White / Green

ROLLEX MEDICAL NZ LTD

PO Box 204159, Highbrook, Auckland 2161, New Zealand
 Phone: +64 9 274 0707 Fax: +64 9 274 0709
 Email: sales@rollexmedical.co.nz www.rollexmedical.co.nz

ROLLEX MEDICAL (AUST) PTY LTD

PO Box 170 | Mt Kuring-gai 2080 | Sydney, NSW | Australia
 Phone: 1300 880 441 | Fax: 1330 880 451
 Email: sales@rollexmedical.com | www.rollexmedical.com.au