



LSFSF39

BENCHTOP LABORATORY FREEZER
SOLID DOOR 39L WHITE



- ✓ **-18°C to -20°C Operating Temperature Range**
Conforms to MHRA guidelines.
- ✓ **Lock Fitted as Standard**
Maintains security and ensures the unit can only be opened by authorised personnel.
- ✓ **Wire Shelves and Basket**
Enables optimum air flow inside the unit.
- ✓ **Unique Anti Microbial Handle**
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **3 Year Extended Warranty**
When serviced by an approved service technician, the warranty is extended to 3 years.



**The Pharmacy
Guild of Australia**



Specifications

Capacity	39L Gross / 35L Net	Spark Free	Yes
Dimensions External	480w x 505d x 565h	Internal Fan	No
Dimensions Internal	360w x 340d x 395h	Temperature range	-18°C – -20°C
Lock	Yes, 2 Keys	Controller	Digital, External
Door	Solid	Interior	HIPS
Door Style	Right Hand Hinge (Reversible)	Housing	Powder coated sheet
Weight	18 kg	Voltage 50Hz	230V
Shelves fit/max	1	Internal light	No
Shelf Material	Wire	Door Open Alarm	Yes
Defrost	Manual	Access Port	Yes
Refrigerant	R600A	Noise Level (db)	-
Adjustable Feet (No & Position)	1 adjustable front foot	Anti Microbial Handle	Yes
Country of Origin	United Kingdom	Colour	White / Green

ROLLEX MEDICAL NZ LTD

PO Box 204159, Highbrook, Auckland 2161, New Zealand
Phone: +64 9 274 0707 Fax: +64 9 274 0709
Email: sales@rollexmedical.co.nz www.rollexmedical.co.nz

ROLLEX MEDICAL (AUST) PTY LTD

PO Box 170 | Mt Kuring-gai 2080 | Sydney, NSW | Australia
Phone: 1300 880 441 | Fax: 1330 880 451
Email: sales@rollexmedical.com | www.rollexmedical.com.au