



# MATOS EXCEL - MEG400

FREESTANDING PHARMACY FRIDGE  
GLASS DOOR 400L WHITE

- ✓ **+2°C to +8°C Operating Temperature Range**  
Conforms to MHRA guidelines.
- ✓ **Lock Fitted as Standard**  
Maintains security and ensures the unit can only be opened by authorised personnel.
- ✓ **Solid or Glass Door Options**  
Most models come with solid or glass door options enabling you to select the format which best suits your requirements.
- ✓ **Wire Shelves**  
Enables optimum air flow inside the unit.
- ✓ **Unique Anti Microbial Handle**  
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **3 Year Extended Warranty**  
When serviced by an approved service technician, the warranty is extended to 3 years.
- ✓ **Min / Max Temperature**  
Contoller has a simple, Read, Record, Reset with the touch of one button.



## Specifications

Capacity	408L Gross / 400L Net	Power	2.6 kWh
Unit Dimensions External	595w x 695d x 1845h	Internal Fan	Yes
Packaged Dimensions	650w x 750d x 1910h	Temperature range	2°C – +8°C
Lock	Yes, 2 Keys	Controller	Digital, External
Door	Glass	Interior	HIPS
Door Style	Right Hand Hinge	Housing	Powder coated sheet
Weight	84 kg	Voltage 50/60Hz	220-240V
Shelves fit/max	6 (5 x full and 1 smaller)	Internal light	Yes
Shelf Material	Wire	Door Open Alarm	Yes
Defrost	Auto	Access Port	Yes - Rear Middle
Refrigerant	R600A	Noise Level (db)	50 dbA
Adjustable Feet (No & Position)	2 adjustable feet	Anti Microbial Handle	Yes
Input Power (W)	240	Colour	White / Green

### ROLLEX MEDICAL NZ LTD

PO Box 204159, Highbrook, Auckland 2161, New Zealand  
Phone: +64 9 274 0707 Fax: +64 9 274 0709  
Email: sales@rollexmedical.co.nz www.rollexmedical.co.nz

### ROLLEX MEDICAL (AUST) PTY LTD

PO Box 170 | Mt Kuring-gai 2080 | Sydney, NSW | Australia  
Phone: 1300 880 441 | Fax: 1330 880 451  
Email: sales@rollexmedical.com | www.rollexmedical.com.au