



# MATOS EXCEL - MES500

## FREESTANDING PHARMACY FRIDGE SOLID DOOR 500L WHITE



- ✓ **+2°C to +8°C Operating Temperature Range**  
Conforms to MHRA guidelines.
- ✓ **Lock Fitted as Standard**  
Maintains security and ensures the unit can only be opened by authorised personnel.
- ✓ **Solid or Glass Door Options**  
Most models come with solid or glass door options enabling you to select the format which best suits your requirements.
- ✓ **Wire Shelves**  
Enables optimum air flow inside the unit.
- ✓ **Unique Anti Microbial Handle**  
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **3 Year Extended Warranty**  
When serviced by an approved service technician, the warranty is extended to 3 years.
- ✓ **Min / Max Temperature**  
Contoller has a simple, Read, Record, Reset with the touch of one button.

### Specifications

Capacity	495L Gross / 480L Net	Power	3.9 kWh
Unit Dimensions External	1180w x 609d x 1530h	Internal Fan	Yes
Packaged Dimensions	1230w x 655d x 1710h	Temperature range	2°C – +8°C
Lock	Yes, 2 Keys	Controller	Digital, External
Door	Solid - Double Door	Interior	HIPS
Door Style	Right Hand Hinge	Housing	Powder coated sheet
Weight	130 kg	Voltage 50Hz	220-240V
Shelves fit/max	8 (6 x full and 2 base)	Internal light	Yes
Shelf Material	Wire	Door Open Alarm	Yes
Defrost	Auto	Access Port	Yes - Rear Middle
Refrigerant	R600A	Noise Level (db)	50 dbA
Lockable Castors	4	Anti Microbial Handle	Yes
Input Power (W)	380	Colour	White / Green

#### ROLLEX MEDICAL NZ LTD

PO Box 204159, Highbrook, Auckland 2161, New Zealand  
Phone: +64 9 274 0707 Fax: +64 9 274 0709  
Email: sales@rollexmedical.co.nz www.rollexmedical.co.nz

#### ROLLEX MEDICAL (AUST) PTY LTD

PO Box 170 | Mt Kuring-gai 2080 | Sydney, NSW | Australia  
Phone: 1300 880 441 | Fax: 1330 880 451  
Email: sales@rollexmedical.com | www.rollexmedical.com.au