

ROLLEX MEDICAL NZ LTD  
 44c Crooks Rd | East Tamaki | Auckland 2013 | New Zealand  
 PO Box 204159 | Highbrook | Auckland 2161 | New Zealand  
 Phone: +64 9 274 0707 | Fax: +64 9 274 0709  
 Email: sales@rollexmedical.co.nz | www.rollexmedical.co.nz

*The Matos® Cloud freezer series is a new range of purpose built Medical Freezers utilising patented technology. The cabinet and componentry have been specially designed for use in the Medical field.*

## Model Variations

- Cloud 300 F – Solid Door
- Cloud 300 F/G – Glass Door



MATOS Cloud 300 F

## Built-In Monitoring capability

Continuous temperature monitoring with sms and email alerts ensures a robust complete monitoring record of your product (See over the page) – even with loss of power.

## Lock Fitted as Standard

Maintains security and ensures the unit can only be opened by authorised personnel.

## Interior Finish

All Matos® Cloud freezers come with a high quality stainless steel interior and use polyurethane foam thermal insulation.

## Solid or Glass Door Options

Most models come with solid or glass door options enabling you to select the format which best suits your requirements.

## Energy Efficient

Matos® Cloud freezers are more energy efficient and kinder to the environment, meaning the running costs of the freezer are kept as low as possible.

## 3 Year Extended Warranty

When serviced by an approved service technician, the warranty is extended to 3 years.

## Specifications

Capacity	300L	Air convection	Natural
Dimensions External	775w x 1735h x 800d	Temperature range	-35°C – 0°C
Dimensions Internal	496w x 1120h x 524d	Controller	Microprocessor with external LCD graphic display
Max shelf workload	10kg	Interior	Stainless Steel
Max unit workload	80kg	Housing	Powder coated sheet (s/steel optional)
Door	Solid/Glass option	Colour	Silver/Grey
Weight	185kg	Temperature accuracy	1°C
Shelves fit/max	3/6	Temperature fluctuation at -10°C	+/- 0.7°C
Web monitoring	Yes	Interface Port	Yes
DAQ battery type (x 2)	12V sealed lead acid (12V 3.5Ah)	LAN Connection	Yes
Refrigerant	R404A	Door Open Alarm	Yes
Refrigerant amount	0.3kg	Access Port	Yes
Castors	No	Internal Light	No
Power	450W	Remote Alarm Contacts	Yes
Voltage 50/60Hz	115/230V	Door Lock	Yes