

ROLLEX MEDICAL NZ LTD

38E Highbrook Drive, East Tamaki, Auckland 2013, New Zealand
PO Box 204159, Highbrook, Auckland 2161, New Zealand
Phone: +64 9 274 0707 Fax: +64 9 274 0709
Email: sales@rollexmedical.co.nz www.rollexmedical.co.nz

The Matos® Plus Eco fridge series is a new range of purpose built Medical and Vaccination Refrigerators utilising patented technology. The cabinet and componentry have been specially designed for use in the Medical field.

Model Variations

Eco 1460 R/DT – Solid Door; hot gas defrost

Eco 1460 R/GDT – Glass Door; hot gas defrost



MATOS Eco 1460 R/GDT

Interior Finish

All Matos® Plus Eco fridges come with a high quality aluminium interior and use polyurethane foam thermal insulation. A stainless steel interior is also offered as an option.

Lock Fitted as Standard

Maintains security and ensures the unit can only be opened by authorised personnel.

Solid or Glass Door Options

Most models come with solid or glass door options enabling you to select the format which best suits your requirements.

Auto Defrost

The advanced controller initializes a regular defrost sequence to decrease the chance of ice build-up while maintaining the temperature below 8°C. Matos® Plus Eco fridges are ordered with a gas defrost capability.

Energy Efficient

Matos® Plus Eco fridges are more energy efficient and kinder to the environment, meaning the running costs of the refrigerator are kept as low as possible.

3 Year Extended Warranty

When serviced by an approved service technician, the warranty is extended to 3 years.

Specifications

| | | | |
|---------------------|----------------------|--------------------------------|--|
| Capacity | 1460L | Power | 550W |
| Dimensions External | 1450w x 1970h x 950d | Air convection | Forced Air |
| Dimensions Internal | 1340w x 1460h x 750d | Temperature range | 0°C – +10°C |
| Max shelf workload | 30kg | Controller | Microprocessor with external LCD graphic display |
| Max unit workload | 300kg | Interior | Aluminium (s/steel optional) |
| Door | Solid/Glass option | Housing | Powder coated sheet (s/steel optional) |
| Weight | 200kg | Temperature resolution | 0.1°C |
| Shelves fit/max | 2 x 3/1 l | Temperature fluctuation at 4°C | +/- 0.5°C |
| Gas defrost | Optional | Voltage 50/60Hz | 115/230V |
| Refrigerant | R404A | Internal light | Yes |
| Refrigerant amount | 0.27kg | Door Open Alarm | Yes |
| Castors | Yes | Access port | Yes |