

## ROLLEX MEDICAL NZ LTD

38E Highbrook Drive, East Tamaki, Auckland 2013, New Zealand  
PO Box 204159, Highbrook, Auckland 2161, New Zealand  
Phone: +64 9 274 0707 Fax: +64 9 274 0709  
Email: sales@rollexmedical.co.nz www.rollexmedical.co.nz

*The Matos® Plus Eco Cooled Incubators are a new range of purpose built Medical Chambers utilising patented technology. The cabinet and componentry have been specially designed for use in the Medical field.*

## Model Variations

Eco 1460 S – Solid Door

Eco 1460 S/G – Glass Door



MATOS Eco 1460 S/G

## Interior Finish

All Matos® Plus Eco **Cooled Incubators** come with a high quality aluminium interior and use polyurethane foam thermal insulation. A stainless steel interior is also offered as an option.

## Lock Fitted as Standard

Maintains security and ensures the unit can only be opened by authorised personnel.

## Solid or Glass Door Options

Most models come with solid or glass door options enabling you to select the format which best suits your requirements.

## Auto Defrost

The advanced controller initializes a regular defrost sequence to decrease the chance of ice build-up while maintaining the temperature.

## Energy Efficient

Matos® Plus Eco **Cooled Incubators** are more energy efficient and kinder to the environment, meaning the running costs of the refrigerator are kept as low as possible.

## 3 Year Extended Warranty

When serviced by an approved service technician, the warranty is extended to 3 years.

## Specifications

Capacity	1460L	Power	550W
Dimensions External	1450w x 1970h x 950d	Air convection	Forced Air
Dimensions Internal	1340w x 1460h x 750d	Temperature range	+3°C – +40°C
Max shelf workload	30kg	Controller	Microprocessor with external LCD graphic display
Max unit workload	300kg	Interior	Aluminium (s/steel optional)
Door	Solid/Glass option	Housing	Powder coated sheet (s/steel optional)
Weight	200kg	Temperature resolution	0.1°C
Shelves fit/max	2 x 3/11	Temperature fluctuation at 4°C	+/- 0.5°C
Gas defrost	Optional	Voltage 50/60Hz	115/230V
Refrigerant	R404A	Door Open Alarm	Yes
Refrigerant amount	0.27kg	Access Port	Yes
Castors	Yes	Internal light	Yes