



## UNDER BENCH PHARMACY FRIDGE SOLID DOOR - 117L

**+2°C to +8°C Operating Temperature Range**

Conforms to MHRA guidelines.

**Lock Fitted as Standard**

Maintains security and ensures the unit can only be opened by authorised personnel.

**Solid or Glass Door Options**

Most models come with solid or glass door options enabling you to select the format which best suits your requirements.

**Wire Shelves**

Enables optimum air flow inside the unit.

**Port Hole Feature as standard**

Access port is available as standard with an internal diameter of 12mm making it easy to fit an external probe for real time monitoring.

**3 Year Extended Warranty**

When serviced by an approved service technician, the warranty is extended to 3 years.

### Specifications

Capacity	117L	Power	0.9 kW/24h
Dimensions External	600w x 600d x 835h	Internal Fan	Yes
Dimensions Internal	480w x 410d x 690h	Temperature range	2°C – +8°C
Lock	Yes, 2 Keys	Controller	Carel - Digital, External
Door	Solid	Interior	HIPS
Door Style	Right Hand Hinge	Housing	Powder coated sheet
Weight	51 kg	Voltage 50/60Hz	230V
Shelves fit/max	2 and 1/2 shelf	Internal light	Yes
Shelf Material	Wire	Door Open Alarm	Yes
Defrost	Auto	Access Port	Yes 12mm
Refrigerant	R600A	Power Watt	107
Feet	4 levelling feet	Insulation (mm)	50
Country of Origin	Denmark	Colour	Silver / White

#### ROLLEX MEDICAL NZ LTD

PO Box 204159, Highbrook, Auckland 2161, New Zealand  
Phone: +64 9 274 0707 Fax: +64 9 274 0709  
Email: sales@rollexmedical.co.nz www.rollexmedical.co.nz

#### ROLLEX MEDICAL (AUST) PTY LTD

PO Box 170 | Mt Kuring-gai 2080 | Sydney, NSW | Australia  
Phone: 1300 880 441 | Fax: 1330 880 451  
Email: sales@rollexmedical.com | www.rollexmedical.com.au